



DEVICE DETAILS

NAME OF DEVICE	3D Medical Face Mask
ESTABLISHMENT NAME	XU CARE TRADING
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA2507223-116544
BRAND NAME	Dr. Chu
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Face masks are made of 3-ply of water resistant outer layer, non woven filter to filter out the dust and germs, excellent filtering melt blown fabric with electrostatic charge, soft and fine inner non woven skin friendly and suitable for sensitive skin. They have bacterial filtration efficiency and particles filtration efficiency for use by humans in any situation such as pandemic, medical sectors and common use for health concern. A new 3D design for easy talking while using the face mask.
DEVICE INTENDED PURPOSE	3 layers filtering face mask respirator and used as to reduce or prevent the transmission of disease. It's single use product and is not recommended for multiple use.
VALIDITY DATE OF REGISTRATION	24/01/2023 - 23/01/2028

LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	3D Medical Face Mask	Adult
2	3D Medical Face Mask	Kids