



DEVICE DETAILS

| | |
|-------------------------------|--|
| NAME OF DEVICE | NEUTROVIS KF94 KOREAN PREMIUM HIJAB FACE RESPIRATOR |
| ESTABLISHMENT NAME | NEUTROVIS SDN BHD |
| ROLE OF ESTABLISHMENT | MANUFACTURER |
| REGISTRATION NO | GA7884323-116592 |
| BRAND NAME | NEUTROVIS |
| MEDICAL DEVICE CATEGORY | SINGLE-USE DEVICES |
| DEVICE GROUPING TYPE | FAMILY |
| DEVICE DESCRIPTION | Neutrovis KF94 Korean Premium Hijab Face Respirator is designed to limit unfiltered air entry with the flaps staying flush against your face and provides a secure seal and double protection from all sides of the face mask. The 2 layers of inner melt-blown fabric and its innovative 3D design allow for more space between your mouth and the face mask and reduced the risk of the wearer readjusting their face mask |
| DEVICE INTENDED PURPOSE | Neutrovis KF94 Korean Premium Hijab Face Respirator with $\geq 98\%$ PFE and $\geq 98\%$ BFE able to filter out bacteria, pollution, pollen and dust particles from our surroundings |
| VALIDITY DATE OF REGISTRATION | 24/01/2023 - 23/01/2028 |

LIST OF DEVICE

| NO | NAME OF DEVICE | IDENTIFIER |
|-------------------|----------------|------------|
| No results found. | | |