



## DEVICE DETAILS

NAME OF DEVICE	HIMAYA 4-PLY PREMIUM HIJAB MEDICAL FACE MASK
ESTABLISHMENT NAME	NEUTROVIS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA7919823-116596
BRAND NAME	HIMAYA
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Himaya 4ply Premium Hijab Medical Face Mask is made from premium medical grade melt-blown filter fabric that will protect you from germs, bacteria and viruses. The soft, non-woven, inner material of the mask provides better comfort and reduced irritation
DEVICE INTENDED PURPOSE	Himaya 4ply Premium Hijab Medical Face Mask with BFE $\geq 98\%$ and PFE $\geq 98\%$ able to filter out live bacteria (viable particles) from 5 to 1 micron and micro particles (nonviable particles) from 1 to 0.1 micron
VALIDITY DATE OF REGISTRATION	24/01/2023 - 23/01/2028

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Himaya 4-Ply Premium Hijab Medical Face Mask	HMY011
2	Himaya 4-Ply Premium Hijab Medical Face Mask	HMY012
3	Himaya 4-Ply Premium Hijab Medical Face Mask	HMY013
4	Himaya 4-Ply Premium Hijab Medical Face Mask	HMY020