




DEVICE DETAILS

NAME OF DEVICE	DISPOSABLE FACE MASK
ESTABLISHMENT NAME	CREATIVE CONTRACT (M) SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GMD69294338517A
BRAND NAME	MEDICOS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	- 3-ply Super Soft PP Non-Woven - Latex-Free & Hypoallergenic - Bacterial Filtration Efficiency > 99% - Particle Filtration Efficiency 0.1%u03BC > 98% - Fluid Resistant - Eye-Shield Visor attachment option is available
DEVICE INTENDED PURPOSE	- Use to cover wearer's mouth and nose against bacteria - Serve as a barrier against contaminants from and to the environment and surroundings - Regarded as a health practice in healthcare to against and prevent the spread of diseases
VALIDITY DATE OF REGISTRATION	29/11/2022 - 28/11/2027

LIST OF DEVICE

NO	NAME OF DEVICE 	IDENTIFIER
1	Fluid Resistant Surgical Face Mask	3066, 3066, 2052, 3052
2	Hijab Sub-micron Surgical Face Mask	105x
3	Rainbow Surgical Face Mask	3050
4	Rainbow Children Surgical Face Mask	3200
5	Children Surgical Face Mask	3202