



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	DISPOSABLE FACE MASK
<b>ESTABLISHMENT NAME</b>	CREATIVE CONTRACT (M) SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GMD69294338517A
<b>BRAND NAME</b>	MEDICOS
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	- 3-ply Super Soft PP Non-Woven - Latex-Free & Hypoallergenic - Bacterial Filtration Efficiency > 99% - Particle Filtration Efficiency 0.1µm > 98% - Fluid Resistant - Eye-Shield Visor attachment option is available
<b>DEVICE INTENDED PURPOSE</b>	- Use to cover wearer's mouth and nose against bacteria - Serve as a barrier against contaminants from and to the environment and surroundings - Regarded as a health practice in healthcare to against and prevent the spread of diseases
<b>VALIDITY DATE OF REGISTRATION</b>	29/11/2022 - 28/11/2027

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Fluid Resistant Surgical Face Mask	3066, 3066, 2052, 3052
2	Hijab Sub-micron Surgical Face Mask	105x
3	Rainbow Surgical Face Mask	3050
4	Rainbow Children Surgical Face Mask	3200
5	Children Surgical Face Mask	3202