



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	FRONTIFIE KF94 FACE RESPIRATOR
<b>ESTABLISHMENT NAME</b>	NEUTROVIS SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA5775223-128738
<b>BRAND NAME</b>	FRONTIFIE
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Frontifie KF94 Face Respirator is designed to limit unfiltered air entry with the flaps staying flush against your face and provides a secure seal and double protection from all sides of the face mask. The 2 layers of inner melt-blown fabric and its innovative 3D design allow for more space between your mouth and the face mask and reduced the risk of the wearer readjusting their face mask
<b>DEVICE INTENDED PURPOSE</b>	Frontifie KF94 Face Respirator with BFE $\geq$ 98% and PFE $\geq$ 98% is 4-Ply Disposable Medical Face mask is to be use as a barrier to prevent air-borne transfer of bacteria and particles. It is designed to be use in hospital, clinic, operating theater and clean room environment , etc
<b>VALIDITY DATE OF REGISTRATION</b>	18/04/2023 - 17/04/2028

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	Frontifie KF94 Face Respirator - Black	FTF005
2	Frontifie KF94 Face Respirator - White	FTF006