



DEVICE DETAILS

NAME OF DEVICE	FRONTIFIE KF94 FACE RESPIRATOR
ESTABLISHMENT NAME	NEUTROVIS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA5775223-128738
BRAND NAME	FRONTIFIE
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Frontifie KF94 Face Respirator is designed to limit unfiltered air entry with the flaps staying flush against your face and provides a secure seal and double protection from all sides of the face mask. The 2 layers of inner melt-blown fabric and its innovative 3D design allow for more space between your mouth and the face mask and reduced the risk of the wearer readjusting their face mask
DEVICE INTENDED PURPOSE	Frontifie KF94 Face Respirator with BFE $\geq 98\%$ and PFE $\geq 98\%$ is 4-Ply Disposable Medical Face mask is to be use as a barrier to prevent air-borne transfer of bacteria and particles. It is designed to be use in hospital, clinic, operating theater and clean room environment , etc
VALIDITY DATE OF REGISTRATION	18/04/2023 - 17/04/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Frontifie KF94 Face Respirator - Black	FTF005
2	Frontifie KF94 Face Respirator - White	FTF006