



DEVICE DETAILS

NAME OF DEVICE	FRONTIFIE 3-PLY PREMIUM MEDICAL FACE MASK
ESTABLISHMENT NAME	NEUTROVIS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA5902723-129108
BRAND NAME	FRONTIFIE
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	This single use medical face mask with BFE \geq 98% and PFE \geq 98% is made from premium medical grade melt-blown filter fabric that will protect you from germs, bacteria and viruses. The soft, non-woven, inner material of the mask provides better comfort and reduced irritation. It is able to filter out live bacteria (viable particles) from 5 to 1 micron and micro particles (nonviable particles) from 1 to 0.1 micron
DEVICE INTENDED PURPOSE	THIS 3 PLY DISPOSABLE MEDICAL AND SURGICAL FACE MASK IS A NON-WOVEN MATERIAL THAT OFFERS FILTRATION-EFFICIENT, FLUID-RESISTANT, SUPERIOR COMFORT AND BREATHABILITY TO HELP REDUCE WEARER TO EXPOSURE OF POTENTIALLY INFECTIOUS PARTICLES, BLOOD AND BODY FLUIDS . IT IS DESIGNED TO BE USED IN HOSPITAL ESPECIALLY IN OPERATING THEATER AND CLEAN ROOM ENVIRONMENT.
VALIDITY DATE OF REGISTRATION	29/04/2023 - 28/04/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Frontifie 3-Ply Premium Medical Face Mask - Blush Pink	OTR003
2	Frontifie 3-Ply Premium Medical Face Mask - Iris Purple	OTR005