



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	FINGERTIP PULSE OXIMETER
<b>ESTABLISHMENT NAME</b>	SYMBIOMED SDN. BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GC2456023-133310
<b>BRAND NAME</b>	ALVITAL
<b>MEDICAL DEVICE CATEGORY</b>	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	The Fingertip Pulse Oximeter is intended for measuring functional oxygen saturation of arterial hemoglobin (SpO <sub>2</sub> ) and pulse rate as non-invasive spot checking in home and professional caring environment. The Fingertip Pulse Oximeter is as a stand-only device for multiple use with 1 AAA Alkaline battery. It includes an OLED display screen and a warning or indicator function for abnormal readings.
<b>DEVICE INTENDED PURPOSE</b>	The Fingertip Pulse Oximeter is intended for measuring functional oxygen saturation of arterial hemoglobin (SpO <sub>2</sub> ) and pulse rate for both adults and adolescent as non-invasive spot checking in home and professional caring environment. It is designed for fingers between 0.8cm and 2.3cm (0,3 inches ~0,9 inches) and for patients during no-motion condition.
<b>VALIDITY DATE OF REGISTRATION</b>	29/05/2023 - 28/05/2028

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	Alvital Fingertip Pulse oximeter	AT101B
2	Alvital Fingertip Pulse oximeter	AT101C
3	Alvital Fingertip Pulse oximeter (Wireless)	AT101I