



## DEVICE DETAILS

NAME OF DEVICE	FINGERTIP PULSE OXIMETER
ESTABLISHMENT NAME	SYMBIOMED SDN. BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC2456023-133310
BRAND NAME	ALVITAL
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The Fingertip Pulse Oximeter is intended for measuring functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate as non-invasive spot checking in home and professional caring environment. The Fingertip Pulse Oximeter is as a stand-only device for multiple use with 1 AAA Alkaline battery. It includes an OLED display screen and a warning or indicator function for abnormal readings.
DEVICE INTENDED PURPOSE	The Fingertip Pulse Oximeter is intended for measuring functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate for both adults and adolescent as non-invasive spot checking in home and professional caring environment. It is designed for fingers between 0.8cm and 2.3cm (0,3 inches ~0,9 inches) and for patients during no-motion condition.
VALIDITY DATE OF REGISTRATION	29/05/2023 - 28/05/2028



## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Alvital Fingertip Pulse oximeter	AT101B
2	Alvital Fingertip Pulse oximeter	AT101C
3	Alvital Fingertip Pulse oximeter (Wireless)	AT101I