



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	SURGICAL FACE MASK
<b>ESTABLISHMENT NAME</b>	JINGGA ANGGUN SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA6305621-78604
<b>BRAND NAME</b>	EMPRO
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	MODEL : VS-913/MM-139/SM-8768/SM-149S/SM-149C/MS31-3P/GS-319/SK-3187 USA ASTM LEVEL 3 (HIGHEST) ASTM F2101-19 BACTERIAL FILTRATION EFFICIENCY: >=99.9% VIRAL FILTRATION EFFICIENCY>=99.9% PARTICLE FILTRATION EFFICIENCY : >= 98% FLUID RESISTANCE:>= 160mmHg BREATHABILITY -DELTA P:<6.0mmH20/cm2 FLAMMABILITY : CLASS 1 PRODUCT IS HYPOALLERGENIC,LATEX AND FIBERGLASS-FREE RECOMMENDED FOR SENSITIVE SKIN SOFT & LIGHTWEIGHT
<b>DEVICE INTENDED PURPOSE</b>	This product is intended for general purpose infection control practices. When properly worn, it reduces potential contact by the wearer to blood and bodily fluids
<b>VALIDITY DATE OF REGISTRATION</b>	21/11/2021 - 20/11/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	VS-913
2	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	SK-3187
3	SURGICAL FACE MASK - M SERIES	MM-139
4	SURGICAL FACE MASK - S SERIES	SM-8768
5	SURGICAL FACE MASK - R SERIES	SM-149C
6	SURGICAL FACE MASK - C SERIES	MS31-3P
7	SURGICAL FACE MASK - G SERIES	GS-319
8	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	TS-8778

NO	NAME OF DEVICE	IDENTIFIER
9	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	PR-1991
10	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	RE3113
11	Surgical Face Mask -S series	SI-8768



**MEDICAL DEVICE AUTHORITY**