



DEVICE DETAILS

NAME OF DEVICE	SURGICAL FACE MASK
ESTABLISHMENT NAME	JINGGA ANGGUN SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA6305621-78604
BRAND NAME	EMPRO
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	MODEL : VS-913/MM-139/SM-8768/SM-149S/SM-149C/MS31-3P/GS-319/SK-3187 USA ASTM LEVEL 3 (HIGHEST) ASTM F2101-19 BACTERIAL FILTRATION EFFICIENCY: >=99.9% VIRAL FILTRATION EFFICIENCY>=99.9% PARTICLE FILTRATION EFFICIENCY : >= 98% FLUID RESISTANCE:>= 160mmHg BREATHABILITY -DELTA P:<6.0mmH20/cm2 FLAMMABILITY : CLASS 1 PRODUCT IS HYPOALLERGENIC,LATEX AND FIBERGLASS-FREE RECOMMENDED FOR SENSITIVE SKIN SOFT & LIGHTWEIGHT
DEVICE INTENDED PURPOSE	This product is intended for general purpose infection control practices. When properly worn, it reduces potential contact by the wearer to blood and bodily fluids
VALIDITY DATE OF REGISTRATION	21/11/2021 - 20/11/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	VS-913
2	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	SK-3187
3	SURGICAL FACE MASK - M SERIES	MM-139
4	SURGICAL FACE MASK - S SERIES	SM-8768
5	SURGICAL FACE MASK - R SERIES	SM-149C
6	SURGICAL FACE MASK - C SERIES	MS31-3P
7	SURGICAL FACE MASK - G SERIES	GS-319
8	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	TS-8778

NO	NAME OF DEVICE	IDENTIFIER
9	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	PR-1991
10	SURGICAL FACE MASK (3PLY) -COPPER OXIDE MASK	RE3113
11	Surgical Face Mask -S series	SI-8768



MEDICAL DEVICE AUTHORITY