



## DEVICE DETAILS

NAME OF DEVICE	NEUTROVIS KF94 KOREAN PREMIUM HIJAB FACE RESPIRATOR
ESTABLISHMENT NAME	NEUTROVIS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA7884323-116592
BRAND NAME	NEUTROVIS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Neutrovis KF94 Korean Premium Hijab Face Respirator is designed to limit unfiltered air entry with the flaps staying flush against your face and provides a secure seal and double protection from all sides of the face mask. The 2 layers of inner melt-blown fabric and its innovative 3D design allow for more space between your mouth and the face mask and reduced the risk of the wearer readjusting their face mask
DEVICE INTENDED PURPOSE	Neutrovis KF94 Korean Premium Hijab Face Respirator with $\geq 98\%$ PFE and $\geq 98\%$ BFE able to filter out bacteria, pollution, pollen and dust particles from our surroundings
VALIDITY DATE OF REGISTRATION	24/01/2023 - 23/01/2028

## LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	Neutrovis KF94 Korean Premium Hijab Face Respirator	NTV189
2	Neutrovis KF94 Korean Premium Hijab Face Respirator	NTV190
3	Neutrovis KF94 Korean Premium Hijab Face Respirator	NTV191
4	Neutrovis KF94 Korean Premium Hijab Face Respirator	NTV192
5	Neutrovis KF94 Korean Premium Hijab Face Respirator	NTV193