



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	NEUTROVIS 3-PLY HIJAB MEDICAL FACE MASK
<b>ESTABLISHMENT NAME</b>	NEUTROVIS SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA3695523-116424
<b>BRAND NAME</b>	NEUTROVIS
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	3 PLY DISPOSABLE MEDICAL FACE MASK MADE OF NON-WOVEN MATERIAL THAT OFFERS FILTRATION-EFFICIENT, FLUID-RESISTANT, SUPERIOR COMFORT AND BREATHABILITY TO HELP REDUCE WEARER TO EXPOSURE OF POTENTIALLY INFECTIOUS PARTICLES, BLOOD AND BODY FLUIDS
<b>DEVICE INTENDED PURPOSE</b>	3-PLY DISPOSABLE MEDICAL FACE MASK IS TO BE USE AS A BARRIER TO PREVENT AIR-BORNE TRANSFER OF BACTERIA AND PARTICLES. IT IS DESIGNED TO BE USE IN HOSPITAL, CLINIC, OPERATING THEATER AND CLEAN ROOM ENVIRONMENT , ETC
<b>VALIDITY DATE OF REGISTRATION</b>	20/01/2023 - 19/01/2028

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Neutrovis 3-Ply Hijab Medical Face Mask	NTV18
2	Neutrovis 3-Ply Hijab Medical Face Mask	NTV27
3	Neutrovis 3-Ply Hijab Medical Face Mask	NTV28
4	Neutrovis 3-Ply Hijab Medical Face Mask	NTV29
5	Neutrovis 3-Ply Hijab Medical Face Mask	NTV95
6	Neutrovis 3-Ply Hijab Medical Face Mask	NTV96