



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	DISPOSABLE MEDICAL FACE MASK
<b>ESTABLISHMENT NAME</b>	ETERNAL MEDICAL SUPPLY SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA9797123-137405
<b>BRAND NAME</b>	AMICA
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Disposable medical mask is composed of Tape, Plastic nose clip, Inner layer Non-woven fabric, Middle Melt-blown Non-woven fabric and outer Non-woven fabric. The melt-blown material acts as the filter that stops microbes from entering or exiting the mask. Most surgical masks feature pleats or folds. Three pleats are used to allow the user to expand the mask such that it covers the area from the nose to the chin. The type is ear loop, where a string-like material is attached to the mask and placed behind the ears.
<b>DEVICE INTENDED PURPOSE</b>	It should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations .
<b>VALIDITY DATE OF REGISTRATION</b>	28/06/2023 - 27/06/2028

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	AMICA DISPOSABLE MEDICAL FACE MASK	MFM17595AD
2	AMICA DISPOSABLE MEDICAL FACE MASK	MFM12585C