



## DEVICE DETAILS

NAME OF DEVICE	DISPOSABLE MEDICAL FACE MASK
ESTABLISHMENT NAME	ETERNAL MEDICAL SUPPLY SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA9797123-137405
BRAND NAME	AMICA
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Disposable medical mask is composed of Tape, Plastic nose clip, Inner layer Non-woven fabric, Middle Melt-blown Non-woven fabric and outer Non-woven fabric. The melt-blown material acts as the filter that stops microbes from entering or exiting the mask. Most surgical masks feature pleats or folds. Three pleats are used to allow the user to expand the mask such that it covers the area from the nose to the chin. The type is ear loop, where a string-like material is attached to the mask and placed behind the ears.
DEVICE INTENDED PURPOSE	It should only be used for patients and other persons to reduce the risk of spread of infections particularly in epidemic or pandemic situations .
VALIDITY DATE OF REGISTRATION	28/06/2023 - 27/06/2028

## LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	AMICA DISPOSABLE MEDICAL FACE MASK	MFM17595AD
2	AMICA DISPOSABLE MEDICAL FACE MASK	MFM12585C