



DEVICE DETAILS

NAME OF DEVICE	SOFT CONTACT LENS
ESTABLISHMENT NAME	SOLAVISH SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB4596723-139096
BRAND NAME	IVY LENSE
MEDICAL DEVICE CATEGORY	MD 0105 - Non-active ophthalmologic devices
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Soft Contact Lens is a coloured contact lens with 57 varieties of colours intended for cosmetic purpose to enhance and change the hue of eyes of the wearer. It is meant for short-term contact and reusable. Soft Contact Lens is indicated for adults for cosmetic effects. Soft Contact Lens comes in package containing preservation solutions of boric acid and sodium chloride salt solutions. Soft Contact Lenses is a sterile product. The lens is designed to be attached to the corneal tear layer and the lens is slightly covered by the eyelid. When blinking, the eyelid slides on the surface of the lens, causing the lens to move slightly. The tears provide the necessary lubrication for the cornea and help to expel foreign matter.
DEVICE INTENDED PURPOSE	Intended for cosmetic purpose to enhanced and change the hue of the eyes of the wearer.
VALIDITY DATE OF REGISTRATION	08/07/2023 - 07/07/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Ivy Lense	ME28
2	Ivy Lense	ME29
3	Ivy Lense	ME30
4	Ivy Lense	ME54
5	Ivy Lense	P01
6	Ivy Lense	P02
7	Ivy Lense	P03
8	Ivy Lense	P05
9	Ivy Lense	FA12-1

NO	NAME OF DEVICE	IDENTIFIER
10	Ivy Lense	FA12-2
11	Ivy Lense	Z6
12	Ivy Lense	Z7



MEDICAL DEVICE AUTHORITY