



DEVICE DETAILS

NAME OF DEVICE	DISPOSABLE SURGICAL MASK (DISPOSABLE FACE MASK)
ESTABLISHMENT NAME	NS MEDIK PHARMA SUPPLIES SDN.BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA4602423-140060
BRAND NAME	EYOUNG
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Disposable Surgical Mask (Disposable face masks) is meant to help block large-particle droplets, splashes, sprays, or splatter that may contain germs (viruses and bacteria), keeping it from reaching your mouth and nose. Supplied as Non-Sterile and Single-Use only. This mask is composed of three parts, the mask body, nose clip and straps. The mask body has three layers, which are thermally bonded by non-woven fabric and filtering material as the main material. Earloop mask, according to the wearing style. Made by PP Non-Woven Material and PE material. EAR-LOOP (Type IIR).
DEVICE INTENDED PURPOSE	Disposable Surgical Mask is intended used for worn by health professionals during surgery and certain health care procedures to catch microorganisms shed in liquid droplets and aerosols from the wearer's mouth and nose. The medical personnel are dressed with them in operation so that the blood, body fluid and exudation of patients, which has potential infection, can be obstructed effectively to prevent the bacteria infection.
VALIDITY DATE OF REGISTRATION	13/07/2023 - 12/07/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	DISPOSABLE SURGICAL MASK (DISPOSABLE FACE MASK)	17.5cm x 9.5cm