



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	DISPOSABLE SURGICAL MASK (DISPOSABLE FACE MASK)
<b>ESTABLISHMENT NAME</b>	NS MEDIK PHARMA SUPPLIES SDN.BHD.
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA4602423-140060
<b>BRAND NAME</b>	EYOUNG
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	SINGLE
<b>DEVICE DESCRIPTION</b>	Disposable Surgical Mask (Disposable face masks) is meant to help block large-particle droplets, splashes, sprays, or splatter that may contain germs (viruses and bacteria), keeping it from reaching your mouth and nose. Supplied as Non-Sterile and Single-Use only. This mask is composed of three parts, the mask body, nose clip and straps. The mask body has three layers, which are thermally bonded by non-woven fabric and filtering material as the main material. Earloop mask, according to the wearing style. Made by PP Non-Woven Material and PE material. EAR-LOOP (Type IIR).
<b>DEVICE INTENDED PURPOSE</b>	Disposable Surgical Mask is intended used for worn by health professionals during surgery and certain health care procedures to catch microorganisms shed in liquid droplets and aerosols from the wearer's mouth and nose. The medical personnel are dressed with them in operation so that the blood, body fluid and exudation of patients, which has potential infection, can be obstructed effectively to prevent the bacteria infection.
<b>VALIDITY DATE OF REGISTRATION</b>	13/07/2023 - 12/07/2028

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	DISPOSABLE SURGICAL MASK (DISPOSABLE FACE MASK)	17.5cm x 9.5cm