



DEVICE DETAILS

NAME OF DEVICE	NEUTROVIS KF94 KOREAN PREMIUM FACE RESPIRATOR
ESTABLISHMENT NAME	NEUTROVIS SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA3758923-116456
BRAND NAME	NEUTROVIS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Neutrovis KF94 Korean Premium Face Respirator is designed to limit unfiltered air entry with the flaps staying flush against your face and provides a secure seal and double protection from all sides of the face mask. The 2 layers of inner melt-blown fabric and its innovative 3D design allow for more space between your mouth and the face mask and reduced the risk of the wearer readjusting their face mask
DEVICE INTENDED PURPOSE	Neutrovis KF94 Korean Premium Face Respirator with $\geq 98\%$ PFE and $\geq 98\%$ BFE able to filter out bacteria, pollution, pollen and dust particles from our surroundings
VALIDITY DATE OF REGISTRATION	21/01/2023 - 20/01/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Neutrovis KF94 Korean Premium Face Respirator	NTV75
2	Neutrovis KF94 Korean Premium Face Respirator	NTV76
3	Neutrovis KF94 Korean Premium Face Respirator	NTV77
4	Neutrovis KF94 Korean Premium Face Respirator	NTV78
5	Neutrovis KF94 Korean Premium Face Respirator	NTV83
6	Neutrovis KF94 Korean Premium Face Respirator	NTV84
7	Neutrovis KF94 Korean Premium Face Respirator	NTV85
8	Neutrovis KF94 Korean Premium Face Respirator	NTV86
9	Neutrovis KF94 Korean Premium Face Respirator	NTV87
10	Neutrovis KF94 Korean Premium Face Respirator	NTV88

NO	NAME OF DEVICE	IDENTIFIER
11	Neutrovis KF94 Korean Premium Face Respirator	NTV89
12	Neutrovis KF94 Korean Premium Face Respirator	NTV90
13	Neutrovis KF94 Korean Premium Face Respirator	NTV91
14	Neutrovis KF94 Korean Premium Face Respirator	NTV92
15	Neutrovis KF94 Korean Premium Face Respirator	NTV122
16	Neutrovis KF94 Korean Premium Face Respirator	NTV125
17	Neutrovis KF94 Korean Premium Face Respirator	NTV126
18	Neutrovis KF94 Korean Premium Face Respirator	NTV127
19	Neutrovis KF94 Korean Premium Face Respirator	NTV128
20	Neutrovis KF94 Korean Premium Face Respirator	NTV129
21	Neutrovis KF94 Korean Premium Face Respirator	NTV130
22	Neutrovis KF94 Korean Premium Face Respirator	NTV133
23	Neutrovis KF94 Korean Premium Face Respirator	NTV165
24	Neutrovis KF94 Korean Premium Face Respirator	NTV166
25	Neutrovis KF94 Korean Premium Face Respirator	NTV167
26	Neutrovis KF94 Korean Premium Face Respirator	NTV168
27	Neutrovis KF94 Korean Premium Face Respirator	NTV169
28	Neutrovis KF94 Korean Premium Face Respirator	NTV170
29	Neutrovis KF94 Korean Premium Face Respirator	NTV171
30	Neutrovis KF94 Korean Premium Face Respirator	NTV172
31	Neutrovis KF94 Korean Premium Face Respirator	NTV173
32	Neutrovis KF94 Korean Premium Face Respirator	NTV174
33	Neutrovis KF94 Korean Premium Face Respirator	NTV184
34	Neutrovis KF94 Korean Premium Face Respirator	NTV185
35	Neutrovis KF94 Korean Premium Face Respirator	NTV186
36	Neutrovis KF94 Korean Premium Face Respirator	NTV187

NO	NAME OF DEVICE	IDENTIFIER
37	Neutrovis KF94 Korean Premium Face Respirator	NTV188
38	Neutrovis KF94 Korean Premium Face Respirator	NTV196
39	Neutrovis KF94 Korean Premium Face Respirator	NTV211
40	Neutrovis KF94 Korean Premium Face Respirator	NTV212
41	Neutrovis KF94 Korean Premium Face Respirator	NTV226
42	Neutrovis KF94 Korean Premium Face Respirator	NTV227
43	Neutrovis KF94 Korean Premium Face Respirator	NTV151
44	Neutrovis KF94 Korean Premium Face Respirator	NTV156
45	Neutrovis KF94 Korean Premium Face Respirator	NTV213
46	Neutrovis KF94 Korean Premium Face Respirator	NTV214
47	Neutrovis KF94 Korean Premium Face Respirator	NTV215
48	Neutrovis KF94 Korean Premium Face Respirator	NTV216
49	Neutrovis KF94 Korean Premium Face Respirator	NTV309
50	Neutrovis KF94 Korean Premium Face Respirator	NTV310
51	Neutrovis KF94 Korean Premium Face Respirator	NTV146
52	Neutrovis KF94 Korean Premium Face Respirator	NTV147
53	Neutrovis KF94 Korean Premium Face Respirator	NTV148
54	Neutrovis KF94 Korean Premium Face Respirator	NTV152
55	Neutrovis KF94 Korean Premium Face Respirator	NTV153
56	Neutrovis KF94 Korean Premium Face Respirator	NTV154
57	Neutrovis KF94 Korean Premium Face Respirator	NTV90
58	Neutrovis KF94 Korean Premium Face Respirator	NTV92
59	Neutrovis KF94 Korean Premium Face Respirator	NTV535
60	Neutrovis KF94 Korean Premium Face Respirator	NTV533