



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	ANTI-VIRAL FACE MASK
<b>ESTABLISHMENT NAME</b>	IDS MEDICAL SYSTEMS (M) SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GMD47677826618A
<b>BRAND NAME</b>	RESPOKARE
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	RespoKare Anti-Viral Facemask is a single use disposable, flat-folded mask that pleated to expand into a concertina-shaped. The device intended use is to have a physical barrier to protect users from viruses, as well as a chemical layer, consist of a hydrophilic plastic coating on the outer layer containing 2% citric acid and a second inner layer treated with metal ions of copper and zinc, to inactivate virus upon contacting the mask.
<b>DEVICE INTENDED PURPOSE</b>	this is a fluid-resistant, single use, disposable face mask that actively kills 99.99% of the tested viruses. It helps reduce exposure to bodily fluids and particulate materials. It breaths well and provides all-day comfort.
<b>VALIDITY DATE OF REGISTRATION</b>	01/10/2023 - 30/09/2028

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Anti-Viral Mask-Adult-24pcs	RK-200-2002A-301
2	Anti-Viral Mask-Adults-1pc	RK-200-2001A-301
3	Anti-Viral Mask-Adults-5pcs	RK-200-2001A-300
4	Anti-Viral Mask-Adults-24pcs	RK-200-2001A-319
5	Anti-Viral Mask-Adults-30pcs	RK-200-2001A-308
6	Anti-Viral Mask-Child-1pc	RK-200-3015A-301
7	Anti-Viral Mask-Child-5pcs	RK-200-3015A-300
8	Anti-Viral Mask-Child-24pcs	RK-200-3015A-319
9	Anti-Viral Mask-Child-10pcs	RK-200-2005A-310