



DEVICE DETAILS

NAME OF DEVICE	Surgical Face Mask 4PLY
ESTABLISHMENT NAME	JINGGA ANGGUN SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA4407021-78614
BRAND NAME	EMPRO
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	USA ASTM LEVEL 3 (HIGHEST) ASTM F2101-19 ASH20-073199-01 2021FM0094R01E BACTERIAL FILTRATION EFFICIENCY: >=99.9% PARTICLE FILTRATION EFFICIENCY : >= 98% FLUID RESISTANCE:>= 160mmHg BREATHABILITY -DELTA P:<6.0mmH2O/cm2 FLAMMABILITY : CLASS 1 PRODUCT IS HYPOALLERGENIC,LATEX AND FIBERGLASS-FREE RECOMMENDED FOR SENSITIVE SKIN SOFT & LIGHTWEIGHT
DEVICE INTENDED PURPOSE	This product is intended for general purpose infection control practices. When properly worn, it reduces potential contact by the wearer to blood and bodily fluids
VALIDITY DATE OF REGISTRATION	21/11/2021 - 20/11/2026

LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL FACE MASK	MS13-4P
2	SURGICAL FACE MASK	4V-414S
3	Surgical Face Mask (4PLY)	w8686
4	SURGICAL FACE MASK (4PLY)	HC-4178