



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	Surgical Face Mask 4PLY
<b>ESTABLISHMENT NAME</b>	JINGGA ANGGUN SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	MANUFACTURER
<b>REGISTRATION NO</b>	GA4407021-78614
<b>BRAND NAME</b>	EMPRO
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	USA ASTM LEVEL 3 (HIGHEST) ASTM F2101-19 ASH20-073199-01 2021FM0094R01E BACTERIAL FILTRATION EFFICIENCY: >=99.9% PARTICLE FILTRATION EFFICIENCY : >= 98% FLUID RESISTANCE:>= 160mmHg BREATHABILITY -DELTA P:<6.0mmH20/cm2 FLAMMABILITY : CLASS 1 PRODUCT IS HYPOALLERGENIC,LATEX AND FIBERGLASS-FREE RECOMMENDED FOR SENSITIVE SKIN SOFT & LIGHTWEIGHT
<b>DEVICE INTENDED PURPOSE</b>	This product is intended for general purpose infection control practices. When properly worn, it reduces potential contact by the wearer to blood and bodily fluids
<b>VALIDITY DATE OF REGISTRATION</b>	21/11/2021 - 20/11/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL FACE MASK	MS13-4P
2	SURGICAL FACE MASK	4V-414S
3	Surgical Face Mask (4PLY)	w8686
4	SURGICAL FACE MASK (4PLY)	HC-4178