



DEVICE DETAILS

NAME OF DEVICE	SURGICAL FACE MASK
ESTABLISHMENT NAME	ICON PACKAGING SDN BHD
ROLE OF ESTABLISHMENT	MANUFACTURER
REGISTRATION NO	GA2852220-49344
BRAND NAME	Icon Protective
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A single use,non sterile and latex free device made of 2 layer of non woven and 1 layer of melt blown. Adjustable to fit faces of various shapes and sizes. The middle layer of face mask (Melt blown) which tested for Bacterial Filtration Efficiency that has sufficient protection against haze particles, dust and harmful droplets. It also high breathability, hypoallergenic, soft and lightweight. This device is appropriate to be used by healthcare professional during surgical procedures to protect both the patient and healthcare personnel.
DEVICE INTENDED PURPOSE	A surgical face mask are intended to be worn by personnel to avoid transfer of bacteria, microorganisms and dust. Face mask may also help reduce exposure of your saliva and respiratory secretions to others.
VALIDITY DATE OF REGISTRATION	09/11/2020 - 08/11/2025

LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	3PLY SURGICAL FACE MASK - EARLOOP (TYPE IIR)	3001
2	3PLY SURGICAL FACE MASK - EARLOOP (ZIPLOCK) (TYPE IIR)	3002
3	3PLY SURGICAL FACE MASK - TIE ON (TYPE IIR)	3003
4	3PLY SURGICAL FACE MASK - HIJAB (TYPE IIR)	3006
5	3Ply Surgical Face Mask	Premium