



DEVICE DETAILS

NAME OF DEVICE	ULTRAFORMER III
ESTABLISHMENT NAME	IDS MEDICAL SYSTEMS (M) SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GB5454723-151679
BRAND NAME	ULTRAFORMER III/CLASYS
MEDICAL DEVICE CATEGORY	MD 1100 - GENERAL ACTIVE MEDICAL DEVICES
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	The ULTRAFORMER III is a non-invasive equipment. It is typically computer-controlled and is capable of producing localized heat within tissues or organs. Energy delivered to the patient is via an externally-mounted transducer. It is a safe and effective device for the treatment of dermatologic diseases. The life cycle of ULTRAFORMER III is 5 years
DEVICE INTENDED PURPOSE	"ULTRAFORMER III (UF3-M300)" is a non-invasive therapeutic device for the treatment of ptosis, axillary hyperhidrosis and correction of skin laxity using the principle of coagulation and necrosis
VALIDITY DATE OF REGISTRATION	05/10/2023 - 04/10/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	ULTRAFORMER III (Main Body)	UF3000001
2	Hand-piece ASSY	UFH000002
3	L4-4.5	UFC344501
4	L7-3.0	UFC373001
5	L7-1.5	UFC371501
6	MF2	UFC3MF203
7	MF6	UFC3MF604
8	MF9	UFC3MF904
9	MF13	UFC3MFD04