



DEVICE DETAILS

NAME OF DEVICE	FINGERTIP PULSE OXIMETER
ESTABLISHMENT NAME	REVO HEALTHCARE SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC2160119-29755
BRAND NAME	CHOICEMMED
MEDICAL DEVICE CATEGORY	MD 1302 - Monitoring devices of vital physiological parameters
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The Fingertip Pulse Oximeter is a portable non invasive device used to measure and record the oxygen saturation of arterial hemoglobin (SpO ₂) and pulse rate of users at any moment. It is a small and compact device which is applicable in the hospital and home use with simple and convenient operation. However, it is reminded that the Fingertip Pulse Oximeter is not suitable for continuous blood oxygen monitoring. The measuring range for the SpO ₂ is 70% to 100% with 1% resolution while the measuring range for the pulse rate is 30 to 250 beats per minute (bpm) with a resolution of 1 bpm.
DEVICE INTENDED PURPOSE	The Fingertip Pulse Oximeter is a handheld non-invasive device intended for spot-checking of oxygen saturation of arterial hemoglobin (SpO ₂) and Pulse Rate of adult, adolescent and child patients in hospitals, hospital type facilities and homecare. It is not for continuous monitoring.
VALIDITY DATE OF REGISTRATION	12/06/2024 - 11/06/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Fingertip Pulse Oximeter MD300C2	MD300C2
2	Fingertip Pulse Oximeter MD300C23	MD300C23