



DEVICE DETAILS

NAME OF DEVICE	JEWETT HYPEREXTENSION CORSET
ESTABLISHMENT NAME	SOUTHERN CRESCENT SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GMD86234929918A
BRAND NAME	ORTHOSOFT
MEDICAL DEVICE CATEGORY	REUSABLE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A Jewett Hyperextension corset is a type of spinal Orthosis that stabilizes the upper thoracic to the lower lumbar spine. It is a dynamic Orthosis which places the spine into hyperextension in an effort to offload the anterior aspect of one or more vertebral bodies. Utilizing three points of pressure, the Jewett is diamond-shaped in appearance with pads at the sternum and lower abdomen. Additionally, there is a padded posterior strap which completes the hyperextension force.
DEVICE INTENDED PURPOSE	The Jewett brace is designed to give support to your thoracic and lumbar spine by preventing twisting and flexion (bending forward). It will also keep your back in the correct position while you are healing from your injury and/or surgery. Indications: in vertebral compression fractures In kyphosis After spinal surgery In case of fracture and trauma in the spine
VALIDITY DATE OF REGISTRATION	20/12/2023 - 19/12/2028

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	ORTHOSOFT JEWETT HYPEREXTENSION CORSET (S)	OS9018A-S
2	ORTHOSOFT JEWETT HYPEREXTENSION CORSET (M)	OS9018A-M
3	ORTHOSOFT JEWETT HYPEREXTENSION CORSET (L)	OS9018A-L
4	ORTHOSOFT JEWETT HYPEREXTENSION CORSET FLEX MODEL(S)	OS9018-S
5	ORTHOSOFT JEWETT HYPEREXTENSION CORSET FLEX MODEL(M)	OS9018-M
6	ORTHOSOFT JEWETT HYPEREXTENSION CORSET FLEX MODEL(L)	OS9018-L