



DEVICE DETAILS

NAME OF DEVICE	Alternating Bubble Mattress with Adjustable Pump System
ESTABLISHMENT NAME	PFLEGE MEDICAL SUPPLIES SDN. BHD.
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA3698819-35036
BRAND NAME	PROFIX
MEDICAL DEVICE CATEGORY	HOSPITAL HARDWARE
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Alternating Bubble Mattress with Adjustable Pump System is designed for the treatment and prevention of pressure sore stage in short term and domiciliary therapies. The PVC mattress is composed of 130 7 cm - high bubble cells, which is particularly comfortable. It is fixed to the bed by extra flaps on both top and bottom side. Besides that, it is easy to install, maintain and use due to the function of comfort range adjustment. The pump can be hung to the end of the bed by using 2 hooks provided.
DEVICE INTENDED PURPOSE	Alternating Bubble Mattress with Adjustable Pump System is intended for the treatment and prevention of pressure sore stage in short term and domiciliary therapies.
VALIDITY DATE OF REGISTRATION	03/10/2024 - 02/10/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Alternating Bubble Mattress with Adjustable Pump System	PROFIX