

DEVICE DETAILS

NAME OF DEVICE	Oral Wound Dressing
ESTABLISHMENT NAME	NUTRIXGOLD (M) SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA3125324-162424
BRAND NAME	Curatick
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	A oral bandage or dressing system comprising a mucous-adhesive Side on a Protection Side, wherein the Attachment Side has a strong adhesion when it is applied and hydrated on mucus. It contains water-soluble polymer. It reacts with moisture (water and saliva) in the mouth within a short time and changes into a gel state. Thereafter, sticky physical properties are generated and attached to the mucous membrane. The Protection Side for protection of the Attachment Side and oral wound consists of water-insoluble polymer. The polymer on the protection side is an insoluble in the water and covers the wound to protect the applicable area from the environment (virus, bacteria, saliva, food, etc.) in the oral cavity. One of the main purpose of the oral bandage is the protection of wound area against the environment (virus, bacteria, saliva, food, etc.) that could be stimulus in the oral cavity.
DEVICE INTENDED PURPOSE	Curatick Oral Wound Dressing is used for protecting oral wounds by forming a physical barrier on the wound.
VALIDITY DATE OF REGISTRATION	14/02/2024 - 13/02/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Oral Wound Dressing	CR12