



DEVICE DETAILS

NAME OF DEVICE	Essential KF94 Face Mask
ESTABLISHMENT NAME	MAXWELL MEDICAL SUPPLIES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA1623422-92532
BRAND NAME	Essential
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Essential KF94 face masks are made of 4-ply of water-resistant outer layer, non-woven filter to filter out microorganisms and dust, excellent filtering melt-blown fabric with electrostatic charge, soft and fine inner non-woven skin-friendly, and suitable for sensitive skin. Suitable for use in situations like a pandemic, medical sectors, and common use for health concerns.
DEVICE INTENDED PURPOSE	To be used by medical personnel during surgical or other medical procedures in order to protect both the patient and the operating personnel from the transfer of microorganisms, body fluid, and particulate material transfer.
VALIDITY DATE OF REGISTRATION	17/05/2022 - 16/05/2027

LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	Essential KF94 Face Mask (Black)	0762983188204
2	Essential KF94 Face Mask (Blue)	0762983188198
3	Essential KF94 Face Mask (Grey)	0762983188181
4	Essential KF94 Face Mask (Pink)	0762983188174