



## DEVICE DETAILS

NAME OF DEVICE	FLEXIBLE KARMAN CANNULA
ESTABLISHMENT NAME	BEN WAY ENTERPRISES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA3490724-162960
BRAND NAME	WOMANCARE GLOBAL
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Sizes 4-6mm directly compatible with the Ipas Single Valve Aspirator. All Flexible Karmen cannulae require adapters for use with Ipas MVA Plus® or Ipas Double Valve Aspirator. Disposable and packed sterile. Softer, more flexible cannula may be better option for retroverted or retroflexed uterine position. All cannulas have either lines or dots space at intervals to assist in gauging depth of insertion/retraction. The Flexible Karman cannulae have series of dots - the first at 6cm from the tip and others at 1cm intervals.
DEVICE INTENDED PURPOSE	Flexible Karman Cannula are intended for uterine aspiration/uterine evacuation in obstetrics and gynecologic patients.
VALIDITY DATE OF REGISTRATION	20/02/2024 - 19/02/2029

## LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	Flexible Karman Cannula, size 4mm	FC4
2	Flexible Karman Cannula, size 5mm	FC5
3	Flexible Karman Cannula, size 6mm	FC6
4	Flexible Karman Cannula, size 7mm	FC7
5	Flexible Karman Cannula, size 8mm	FC8
6	Flexible Karman Cannula, size 9mm	FC9
7	Flexible Karman Cannula, size 10mm	FC10
8	Flexible Karman Cannula, size 12mm	FC12