



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	FLEXIBLE KARMAN CANNULA
<b>ESTABLISHMENT NAME</b>	BEN WAY ENTERPRISES SDN BHD
<b>ROLE OF ESTABLISHMENT</b>	AUTHORIZED REPRESENTATIVE
<b>REGISTRATION NO</b>	GA3490724-162960
<b>BRAND NAME</b>	WOMANCARE GLOBAL
<b>MEDICAL DEVICE CATEGORY</b>	SINGLE-USE DEVICES
<b>DEVICE GROUPING TYPE</b>	FAMILY
<b>DEVICE DESCRIPTION</b>	Sizes 4-6mm directly compatible with the Ipas Single Valve Aspirator. All Flexible Karman cannulae require adapters for use with Ipas MVA Plus® or Ipas Double Valve Aspirator. Disposable and packed sterile Softer, more flexible cannula may be better option for retroverted or retroflexed uterine position. All cannulas have either lines or dots space at intervals to assist in gauging depth of insertion/retraction. The Flexible Karman cannulae have series of dots - the first at 6cm from the tip and others at 1cm intervals.
<b>DEVICE INTENDED PURPOSE</b>	Flexible Karman Cannula are intended for uterine aspiration/uterine evacuation in obstetrics and gynecologic patients.
<b>VALIDITY DATE OF REGISTRATION</b>	20/02/2024 - 19/02/2029

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	Flexible Karman Cannula, size 4mm	FC4
2	Flexible Karman Cannula, size 5mm	FC5
3	Flexible Karman Cannula, size 6mm	FC6
4	Flexible Karman Cannula, size 7mm	FC7
5	Flexible Karman Cannula, size 8mm	FC8
6	Flexible Karman Cannula, size 9mm	FC9
7	Flexible Karman Cannula, size 10mm	FC10
8	Flexible Karman Cannula, size 12mm	FC12