



DEVICE DETAILS

NAME OF DEVICE	3-PLY SURGICAL FACE MASK
ESTABLISHMENT NAME	WATSON'S PERSONAL CARE STORES SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GA10955724-172780
BRAND NAME	WATSONS
MEDICAL DEVICE CATEGORY	SINGLE-USE DEVICES
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	A 3-Ply Surgical Face Mask (ASTM Level 3, EN14683 Type IIR) made with 3-ply pleated PP non woven material. Supplied in packs of 10's.
DEVICE INTENDED PURPOSE	Disposable device that can filter droplets, pollens, dust and particulates.
VALIDITY DATE OF REGISTRATION	29/05/2024 - 28/05/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	WATSONS 3-PLY SURGICAL FACE MASK ADULTS 10'S	37506
2	WATSONS 3-PLY SURGICAL FACE MASK KIDS 10'S	37632