



DEVICE DETAILS

NAME OF DEVICE	PULSE OXIMETER
ESTABLISHMENT NAME	JITRON HEALTHCARE SDN BHD
ROLE OF ESTABLISHMENT	AUTHORIZED REPRESENTATIVE
REGISTRATION NO	GC8919924-173487
BRAND NAME	JUMPER
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The fingertip pulse Oximeter features small size, low power consumption, convenient operation, and portability. It is only necessary for a user to put one of his fingers into the fingertip photoelectric sensor for diagnosis, and a display screen will show oxygen saturation. IT IS NOT USED FOR CONTINUOUS MONITORING AND IS NOT USED TOGETHER WITH OTHER MONITORING DEVICES.
DEVICE INTENDED PURPOSE	JPD-500A The Fingertip Pulse Oximeter is non-invasive device intended for spot-checking of functional oxygen saturation of arterial hemoglobin (Spo2) and pulse rate. The portable fingertip device is indicated for adult and pediatric patients in home and hospital environments. JPD-500E The Fingertip Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO2) and the pulse rate of adult and pediatric patients in home and hospital environments.
VALIDITY DATE OF REGISTRATION	05/06/2024 - 04/06/2029

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	JUMPER FINGERTIP PULSE OXIMETER (500E)	JPD-500E
2	JUMPER PULSE OXIMETER JPD-500A	JPD-500A