



DEVICE DETAILS

NAME OF DEVICE	SAMSUNG LABGEO HC10 HEMATOLOGY ANALYZER
ESTABLISHMENT NAME	SILVER STATE LABORATORY SDN BHD
REGISTRATION NO	IVDB10668419-35165
BRAND NAME	SAMSUNG
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	The Samsung LABGEO HC10 Hematology Analyzer measures impedance by using the Coulter Principle for counting cells passing through an aperture and optically measures the haemoglobin content of red blood cells. The analyzer features a graphical color LCD module with a touch screen and a START button. The software interfaces with a 58 mm built-in thermal printer module and allows sending results to an external printer via a USB port.
DEVICE INTENDED PURPOSE	The Samsung LABGEO HC10 is a quantitative, automated haematology analyzer and leukocyte differential counter for in vitro diagnostic use of whole blood specimens, developed for clinical laboratories.
VALIDITY DATE OF REGISTRATION	04/10/2019 - 03/10/2024

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
No results found.		