



DEVICE DETAILS

NAME OF DEVICE	CLEARFLEX ECONOMY SILICONE ANAESTHETIC FACE MASKS
ESTABLISHMENT NAME	KL MED SUPPLIES (M) SDN BHD
REGISTRATION NO	GB441281327219
BRAND NAME	INTERSURGICAL
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	ECONOMY SILICONE ANAESTHETIC MASKS.LATEX FREE SILICONE MASK. AS A FACEMASK TO DELIVER ANAESTHETIC AND RESPIRATORY GASES OR A BRONCHOSCOPY MASK TO DELIVER ANALGESIC GASES.
DEVICE INTENDED PURPOSE	TO DELIVER RESPIRATORY GASES TO A PATIENT VIA A MASK.
VALIDITY DATE OF REGISTRATION	08/04/2019 - 07/04/2024

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	CLEARFLEX ANAESTHESIA MASK, SILICONE, FLEXIBLE FLANGE, SIZE 0	8848000
2	CLEARFLEX ANAESTHESIA MASK, SILICONE, FLEXIBLE FLANGE, SIZE 1	8848001
3	CLEARFLEX ANAESTHESIA MASK, SILICONE, FLEXIBLE FLANGE, SIZE 2	8848002
4	CLEARFLEX ANAESTHESIA MASK, SILICONE, FLEXIBLE FLANGE, SIZE 3	8848003
5	CLEARFLEX ANAESTHESIA MASK, SILICONE, FLEXIBLE FLANGE, SIZE 4	8848004
6	CLEARFLEX ANAESTHESIA MASK, SILICONE, FLEXIBLE FLANGE, SIZE 5	8848005