



## DEVICE DETAILS

|                               |  |
|-------------------------------|--|
| NAME OF DEVICE                | NONIN ONYX VANTAGE FINGER PULSE OXIMETER   |
| ESTABLISHMENT NAME            | ACUCARE SYSTEMS (M) SDN BHD  |
| REGISTRATION NO               | GC11120720-38475   |
| BRAND NAME                    | NONIN  |
| DEVICE GROUPING TYPE          | FAMILY   |
| DEVICE DESCRIPTION            | Noinin Medical's finfer pulse oximeters are small lightweight portable devices that display numerical values for functional oxygen saturation of arterial hemoglobin (SpO2) and pulse by measuring the absorption of red and infrared (IR) light passing through perfused tissue. Changes in the absorption caused by the pulsation of blood in the vascular bed are used to determine oxygen saturation and pulse rate. |
| DEVICE INTENDED PURPOSE       | The Nonin Finger Pulse Oximeter is a small, lightweight portable device intended for used in measuring and displaying functional hemoglobin (%SpO2) and pulse rate of patients who are well or poorly perfused.  |
| VALIDITY DATE OF REGISTRATION | 06/01/2020 - 05/01/2025  |

## LIST OF DEVICE



MEDICAL DEVICE AUTHORITY

| NO | NAME OF DEVICE                   | IDENTIFIER |
|----|----------------------------------|------------|
| 1  | Onyx II 9550 Pulse Oximeter      | 9550       |
| 2  | Onyx Vantage 9590 Pulse Oximeter | 9590       |
| 3  | Onyx3 9591 Finger Pulse Oximeter | 9591       |
| 4  | 9500CC- Carrying Case, Black     | 2003-102   |
| 5  | 9590 Lanyard                     | 6918-000   |
| 6  | Retractable Lanyard              | 6919-002   |
| 7  | Fingertip Carrying Case          | 5097-000   |
| 8  | Protective Carrying Case         | 3802-000   |
| 9  | Fingertip Lanyard                | 5096-000   |