



## DEVICE DETAILS

<b>NAME OF DEVICE</b>	FINGERTIP PULSE OXIMETER
<b>ESTABLISHMENT NAME</b>	GASWORLD SDN BHD
<b>REGISTRATION NO</b>	GC93820154517
<b>DESCRIPTION</b>	PULSE OXIMETER IS A VERY IMPORTANT AND COMMON DEVICE TO CHECK PATIENT'S BLOOD OXYGEN SATURATION (SPO2) LEVEL AND PULSE RATE. AS A PORTABLE, COMPACT, SIMPLE, RELIABLE AND DURABLE PHYSIOLOGICAL MONITORING DEVICE, HANDHELD PULSE OXIMETER AH- MX CAN GREATLY ENHANCES PATIENT CARE. IT IS WIDELY APPLIED IN CLINIC, HOSPITAL, SOCIAL MEDICAL ORGANIZATIONS, AND FIRST AID AND SO ON.
<b>DEVICE GROUPING TYPE</b>	Family
<b>INTENDED PURPOSE</b>	THE AX-XX SERIES PULSE OXIMETER IS INTENDED FOR MEASURING AND DISPLAYING FUNCTIONAL OXYGEN SATURATION OF ARTERIAL HEMOGLOBIN (SPO2) AND PULSE RATE FOR ADULT, PEDIATRIC IN FAMILY, HOSPITAL (INCLUDING CLINICAL USE IN INTERNIST/SURGERY, ANESTHESIA, INTENSIVE CARE AND ETC.), OXYGEN BAR, SOCIAL MEDICAL ORGANIZATIONS, PHYSICAL CARE IN SPORTS(IT CAN BE USED BEFORE OR AFTER SPORTS. CONTINUOUSLY MONITORING IS NOT RECOMMENDED.)
<b>VALIDITY DATE OF REGISTRATION</b>	07-02-2017 - 06-02-2022

## LIST OF DEVICE

#	NAME OF DEVICE	IDENTIFIER
1	Fingertip Pulse Oximeter	AE-O2
2	Pediatric Fingertip Pulse Oximeter (Black Panda)	AE-C9
3	Pediatric Fingertip Pulse Oximeter (Blue Panda)	AE-CA
4	Pediatric Fingertip Pulse Oximeter (Pink Panda)	AE-CB
5	Fingertip Pulse Oximeter	AE-M1