



## DEVICE DETAILS

NAME OF DEVICE	Non-woven face mask
ESTABLISHMENT NAME	BIOD MEDICA SDN. BHD.
REGISTRATION NO	GA8910620-41613
BRAND NAME	HelloDr Face Mask
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Single use 3 ply face mask is made of polypropylene spunbond and meltblown nonwoven fabrics with adjustable nose strip and flat ties or round elastic ear- loops to give basic protection to the user's nose and mouth area.
DEVICE INTENDED PURPOSE	Face mask generally use by but not limited to operating room personnels and other general healthcare workers to provide basic protection to both patients and healthcare workers against transfer of microorganisms, blood and body fluids, and airborne particulates.
VALIDITY DATE OF REGISTRATION	07/04/2020 - 06/04/2025

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Non-woven face mask	MF05A
2	Non-woven face mask	MF05B