



## DEVICE DETAILS

NAME OF DEVICE	FRONTIFIE PREMIUM MEDICAL FACE MASK
ESTABLISHMENT NAME	ROCKY SCIENTIFIC INSTRUMENTATION
REGISTRATION NO	GA2211920-50702
BRAND NAME	FRONTIFIE
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	THE PRODUCT IS MADE OF MELT-BLOWN FABRIC WITH BFE (BACTERIAL FILTRATION EFFICIENCY) >98% AND PFE ( PARTICLE FILTRATION EFFICIENCY) >96%, WHICH IS USED TO BLOCK LIVE BACTERIA (VIABLE PARTICLES) FROM 1 TO 5 MICRON AND MICRO PARTICLES (NONVIABLE PARTICLES) FROM 0.1 MICRON TO 1 MICRON
DEVICE INTENDED PURPOSE	to provide a barrier and to protect the patient and hospital staff or doctors from direct transmission of infective agents, blood,splash,viruses, suitable to be used in a medical, surgical,dental and public
VALIDITY DATE OF REGISTRATION	17/12/2020 - 16/12/2025

## LIST OF DEVICE



MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	FRONTIFIE PREMIUM MEDICAL FACE MASK	OTR003
2	FRONTIFIE PREMIUM MEDICAL FACE MASK	OTR005