



DEVICE DETAILS

NAME OF DEVICE	SURGICAL NON-STERILE DISPOSABLE FACE MASK
ESTABLISHMENT NAME	EUGIN RESOURCES SDN BHD
REGISTRATION NO	GA3930621-60398
BRAND NAME	EUGINMask
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	Surgical mask is a single use 3 layers mask with outer and inner layer (spunbond polypropylene) that sandwich a melt-blown polypropylene filter material. Earloops are of urethane elastic fiber and latex-free. The nose bridge is a pliable white aluminium stripe, covered by PP covering. All of the materials used in the construction of the Surgical Mask are being used in current marketed device.
DEVICE INTENDED PURPOSE	Surgical mask is intended for single use by operating room personnel and other general healthcare workers to protect both patients and healthcare workers against transfer of microorganisms, blood and body fluids, and airborne particulates.
VALIDITY DATE OF REGISTRATION	30/04/2021 - 29/04/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	SURGICAL NON-STERILE DISPOSABLE FACE MASK	SNSDFM