



## DEVICE DETAILS

NAME OF DEVICE	DISPOSABLE MEDICAL FACE MASK
ESTABLISHMENT NAME	SIMPLY K SDN BHD
REGISTRATION NO	GA8323421-65220
BRAND NAME	MEDICMAD
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	Disposable medical face masks comply with EN 14683 'Type IIR' embedded with nose clip and ear loops or head loops. This medical face mask made of 3-ply of premium non-woven fabric with Bacterial Filtration Efficiency (BFE) >98 and Splash Resistance Pressure > 16.0 kPa.
DEVICE INTENDED PURPOSE	Disposable Medical face mask that covers the user's nose and mouth and provides a physical barrier from transmission of viruses, microorganisms and dust. It is intended to be use by healthcare provider and other persons for medical purposes and to control infectious diseases, particularly in epidemic or pandemic situations.
VALIDITY DATE OF REGISTRATION	28/06/2021 - 27/06/2026

## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	DISPOSABLE MEDICAL FACE MASK	ADULT MEDICAL FACE MASK - EAR LOOP
2	DISPOSABLE MEDICAL FACE MASK	ADULT MEDICAL FACE MASK - HEAD LOOP
3	DISPOSABLE MEDICAL FACE MASK	CHILDREN MEDICAL FACE MASK - EAR LOOP
4	DISPOSABLE MEDICAL FACE MASK	CHILDREN MEDICAL FACE MASK - HEAD LOOP