



DEVICE DETAILS

NAME OF DEVICE	Durio 3 Ply Surgical Face Mask
ESTABLISHMENT NAME	DURIO PPE SDN BHD
REGISTRATION NO	GA9895121-66714
DESCRIPTION	Surgical mask is fluid-resistant, disposable, and loose-fitting devices that create a physical barrier between the mouth and nose of the wearer and the immediate environment.
DEVICE GROUPING TYPE	Family
INTENDED PURPOSE	Disposable device that covers the user's nose and mouth and provides a physical barrier to droplets, fluids, bacteria and particulate materials.
VALIDITY DATE OF REGISTRATION	16-07-2021 - 15-07-2026

LIST OF DEVICE

#	NAME OF DEVICE	IDENTIFIER
1	Hijab 3 Ply Surgical Face Mask	521
2	X' Mas 3 Ply Surgical Face Mask	533
3	Raya 3 Ply Surgical Face Mask	534
4	Gong Xi Fa Chai 3 Ply Surgical Face Mask	535