



## DEVICE DETAILS

NAME OF DEVICE	PULSE OXIMETER
ESTABLISHMENT NAME	QUANTUM UPSTREAM SDN.BHD.
REGISTRATION NO	GC10889321-68404
BRAND NAME	EZYSWISS
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	ezyswiss Pulse Oximeter is intended for non-invasive spot-check measurement of functional oxygen saturation of arterial hemoglobin (SpO2). Advanced DSP algorithm can minimize the influence of motion artifact and improve measurement accuracy of low perfusion. The Pulse Oximeter can be used to measure human SpO2 and heart rate through finger. The product is suitable for family, hospital (including clinical use in internist/surgery, Anesthesia, pediatrics etc.) Oxygen Bar, social medical organizations, physical care in sports etc.
DEVICE INTENDED PURPOSE	ezyswiss Pulse Oximeter Pulse Oximeter is a non-invasive device intended for spot checking of functional oxygen saturation of arterial hemoglobin (SpO2) and pulse rate (PR). This portable device is indicated for use in adult patients in clinical institution and home environments.
VALIDITY DATE OF REGISTRATION	03/08/2021 - 02/08/2026



## LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	PULSE OXIMETER	OXI-01