



DEVICE DETAILS

NAME OF DEVICE	WALKING STICK QUAD CANE
ESTABLISHMENT NAME	TEEPHAM MEDICAL SDN. BHD.
REGISTRATION NO	GMD5503922416A
BRAND NAME	NURCARE
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	ALK524 & ALK521 is made of high-strength aluminum alloy tube & four steel ferrules at base. Plastic hand grip is durable and comfortable. Adjustable height with 4 non-slip rubber foot, stable and durable. The quad canes comes in narrow or wide base.
DEVICE INTENDED PURPOSE	Walking stick quad cane is a device used by elderly or disabled people as a walking aid. It has four ferrules at the base, enabling them to stand freely, and offering a more firm base for standing.
VALIDITY DATE OF REGISTRATION	04/10/2021 - 03/10/2026

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Walking stick Quad Cane	ALK524
2	Walking stick Quad Cane	ALK521