



## DEVICE DETAILS

NAME OF DEVICE	ISOLATOR PLUS FACE MASK
ESTABLISHMENT NAME	CANTEL MEDICAL (MALAYSIA) SDN. BHD.
REGISTRATION NO	GMD9353389617A
BRAND NAME	ISOLATOR PLUS
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	The N95 respirator masks consists of a polypropylene meltblown filter media sandwiched between non-woven wetlaid outer and inner covers. The mask is secured using two polyester covered non-latex rubber headbands.
DEVICE INTENDED PURPOSE	The N95 surgical mask is intended to be worn by healthcare personnel during surgical procedures to protect both the pateint and healthcare personnel from the transfer of microorganisms, body fluids, and particulate materials. Intended use includes use as a procedure mask, isolation mask and/or dental face mask.
VALIDITY DATE OF REGISTRATION	16/01/2022 - 15/01/2027

## LIST OF DEVICE



NO	NAME OF DEVICE	IDENTIFIER
1	ISOLATOR PLUS FACE MASK	GPRN95