



DEVICE DETAILS

NAME OF DEVICE	HANDHELD HANDHELD PULSE OXIMETER
ESTABLISHMENT NAME	IHB HEALTHCARE SDN BHD
REGISTRATION NO	GC5337721-78806
BRAND NAME	AEON TECHNOLOGY
DEVICE GROUPING TYPE	SYSTEM
DEVICE DESCRIPTION	Aeon Technology A36 Handheld Handheld Pulse Oximeter is a non-invasive, handheld patient Oximeter. It operates on alkaline or rechargeable battery power supply. It is compact, small, light, and easy for learning and handling. It is suitable for monitoring adult and child patients. It is widely used only in Professional medical organizations, such as, hospitals, clinics. Parameters measured by the oximeter include: arterial oxygen saturation (SpO2), pulse rate (PR), bargraph and plethysmogram. The oximeter measures these parameters through a SpO2 sensor and displays them on the color TFT LCD screen after certain further processing. The oximeter is operated and controlled by the buttons on the front panel. It adopts a 2.8 inch color TFT LCD screen in displaying measurements and in supplementary status indication.
DEVICE INTENDED PURPOSE	Aeon Technology Handheld Pulse Oximeter is a non-invasive device intended for monitoring adult and child patients. It is widely used only in Professional medical organizations, such as, hospitals, clinics.
VALIDITY DATE OF REGISTRATION	21/11/2021 - 20/11/2026

MEDICAL DEVICE AUTHORITY

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	Handheld Pulse Oximeter	A36
2	Adult finger SpO2 sensor	N.A.
3	USB line	N.A.
4	1.5V Alkaline Battery	N.A.
5	Optical disk	N.A.
6	Neonatal finger SpO2 sensor	N.A.
7	Child finger SpO2 sensor	N.A.
8	Temperature probe	N.A.
9	Rechargeable Li battery	N.A.