



## DEVICE DETAILS

NAME OF DEVICE	FLEXOR CHECK-FLO INTRODUCER
ESTABLISHMENT NAME	COOK ASIA (MALAYSIA) SDN BHD
REGISTRATION NO	GB5300419-34174
BRAND NAME	COOK
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	FLEXOR CHECK-FLO INTRODUCERS ARE DESIGNED TO PERFORM AS A GUIDING SHEATH AND/OR INTRODUCER SHEATH, AND CAN BE ORDERED WITH VARIOUS PROXIMAL FITTINGS SUCH AS A TUOHY-BORST VALVE OR A CHECK-FLO® VALVE. SOME INTRODUCER SETS ARE AVAILABLE WITH A RADIOPAQUE BAND INCORPORATED WITHIN THE SHEATH MATERIAL TO IDENTIFY THE LOCATION OF THE SHEATH'S DISTAL TIP, WHILE OTHER SETS ARE AVAILABLE WITH AQ® HYDROPHILIC COATING TO PROVIDE A LUBRICIOUS SURFACE TO EASE INSERTION.
DEVICE INTENDED PURPOSE	FLEXOR CHECK-FLO INTRODUCERS ARE INTENDED FOR INTRODUCTION OF BALLOONS, CLOSED AND NON-TAPERED END CATHETERS OR OTHER DIAGNOSTIC AND INTERVENTIONAL DEVICES.
VALIDITY DATE OF REGISTRATION	29/08/2019 - 28/08/2024

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
No results found.		