



## DEVICE DETAILS

NAME OF DEVICE	NOVOTEMP MULTI-FUNCTION INFRARED THERMOMETER
ESTABLISHMENT NAME	REVO HEALTHCARE SDN BHD
REGISTRATION NO	GB6128319-26799
BRAND NAME	NOVOTEMP
DEVICE GROUPING TYPE	FAMILY
DEVICE DESCRIPTION	The TH50Z and TH520 Infrared Thermometer allows measurement of body temperature. The temperature is taken between the forehead and the temple or alternatively in the ear. The reading is displayed within a few seconds. Apart from measuring body temperature, the TH50Z and TH520 Infrared Thermometer can also be used to measure the temperature of objects (e.g. fluids, surfaces, ambient temperature). The temperature range for this function is between -22°C and 80°C. The reading recorded is displayed within seconds.
DEVICE INTENDED PURPOSE	The TH50Z and TH520 Infrared Thermometer is intended for the intermittent measurement and monitoring of human body temperature by carrying out measurement in the ear or on the central forehead. Additionally, the thermometer is suitable to measure the temperature of objects in the range of -22°C to 80°C (-7.6°F - 176°F).
VALIDITY DATE OF REGISTRATION	13/03/2019 - 12/03/2024

## LIST OF DEVICE

MEDICAL DEVICE AUTHORITY

NO	NAME OF DEVICE	IDENTIFIER
1	TH50Z Multi-function Infrared Thermometer	TH50Z
2	TH520 Multi-function Infrared Thermometer	TH520