



DEVICE DETAILS

NAME OF DEVICE	PULSE OXIMETER CMS50DL2
ESTABLISHMENT NAME	SOLID INTELLECTUAL SDN. BHD.
REGISTRATION NO	GC1687922-85802
BRAND NAME	CMS50DL2
DEVICE GROUPING TYPE	SINGLE
DEVICE DESCRIPTION	The Fingertip Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO2) and the pulse rate of adult and pediatric patients in home and hospital environments (including clinical use in internist/surgery, anesthesia, intensive care etc.). This device is not intended for continuous monitoring. The Pulse Oximeter features in small volume, low power consumption, convenient operation and being portable. It is only necessary for patient to put one of his fingers into a fingertip photoelectric sensor for diagnosis, and a display screen will directly show measured value of Hemoglobin Saturation. The device doesn't incorporate tissues, or blood products. The device is reusable, noninvasive, implantable device.
DEVICE INTENDED PURPOSE	The Fingertip Pulse Oximeter is a non-invasive device intended for the spot-check of oxygen saturation of arterial hemoglobin (SpO2) and the pulse rate of adult and pediatric patients in home and hospital environments (including clinical use in internist/surgery, anesthesia, intensive care ect.). This device is not intended for continuous monitoring.
VALIDITY DATE OF REGISTRATION	27/01/2022 - 26/01/2027

MEDICAL DEVICE AUTHORITY

LIST OF DEVICE

NO	NAME OF DEVICE	IDENTIFIER
1	PULSE OXIMETER CMS50DL2	CMS50DL2